

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922) 49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58
Иваново (4932)77-34-06
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Киргизия (996)312-96-26-47

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Россия (495)268-04-70

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Сургут (3462)77-98-35
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Казахстан (772)734-952-31

Тольяти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

<https://grant.nt-rt.ru> || gtq@nt-rt.ru

Устройства криоконсервации



Cryopreservation

CRF-1

Liquid nitrogen and cryogen free controlled rate freezer for cryopreservation

CRF-1

Controlled Rate Freezer

The CRF-1 liquid nitrogen and cryogen free controlled rate freezer brings accuracy, precision and reproducibility to biological cryopreservation. This unit is ideal for research into the cryopreservation of a wide range of material including embryos, stem cells, mammalian cells, spermatozoa, antibodies, tissue sections and rodent organs.

The CRF-1 has been developed to minimise and reduce the risk of contamination and designed for use in cleanrooms and barrier facilities. The CRF-1 fits neatly on a bench-top. Its performance of cell viability after freezing is equal if not more advanced to comparable liquid nitrogen freezers. As alcohol is not used, there is also a greatly reduced fire risk. The CRF-1 will cool down to -100°C with straws.

The cooling rate of the CRF-1 is precisely controlled, ensuring accuracy and reproducibility throughout the freezing profile, especially for the important nucleation/seeding phase. This provides optimal recovery of cells on thawing. Operation is simple and can be carried out with or without a PC; data can be logged via PC software and cooling profiles are directly displayed on PC screen. Different cooling profiles are available from a drop-down menu and customised profiles can be written.



CRF-1

Controlled Rate Freezer

Main applications

The CRF-1 is highly versatile and can be used for the cryopreservation research of a wide range of samples in cryovials, straws, bags, microplates and Matrix-96-well block plates in the following areas:

- Transgenic embryos research
- Stem cell research
- Clinical and research samples, e.g. lymphocytes and tissue cell lines in conventional cryovials
- Various mammalian cells including cardiomyocytes, adipose, liver and muscle
- Cord blood derived stem cells
- Adherent cells and stem cells in microplates
- Cell suspensions in numbered/barcoded arrays
- Robotic integration – the CRF-1 has also successfully been integrated into robotic systems
- Suitable for applications in veterinarian IVF

Key benefits and features

- Accurate and reproducible control of cooling rates and sample temperatures
- Easy to use and samples can be nucleated/seeded in-situ
- Linear and non-linear cooling profiles
- Low running costs: estimated at 1% of liquid nitrogen-controlled rate freezing
- Temperature remains at -100°C at the end of cycle for straw applications until freezer is switched off
- Un-interruptible Power Supply (UPS): complete cycle run if power fails (supplied as an optional accessory)
- CE marked (laboratory use)
- Servicing and calibration available
- 3 year warranty



CRF-1

Technical specifications



		CRF-1
		Controlled rate freezer
Dimensions	w x d x h mm	264 x 367 x 445
Weight	Kg	9.8
Max plate temperature	°C	30
Min plate temperature	°C	-80 or -100 (dependent on model)
Temperature accuracy	±°C	0.5
Temperature stability	±°C	0.1
Temperature uniformity	±°C	1.0
Max cooling rate		10°C /min
Max temperature stabilisation time		5 minutes

CRF-1

Optional accessories and services









CRF-1 API	API package which unlocks the CRF-1 for use with users own application
CAL CRF-1	Calibration certificate
CRF-1 MAINT 1	1 year maintenance agreement includes clean and test report (excludes repair and replacement of faulty parts)
IQOQ CRF-1	IQOQ Booklet for CRF-1
PQ CRF-1	PQ Booklet for CRF-1
CRF-1 CRYOPEN	Cryopen nucleating tool
CRF-1 CARTRIDGE	Replacement cartridge for Cryopen
CRF-1 UPS	Un-interruptible Power Supply for the CRF-1 for a 3 hour cycle in the event of an electrical power failure. 220-240v operation.

CRF-1

Controlled Rate Freezer

Product range

The range includes various models which provide optimum performance for a specific common vessel or can accommodate two different types vessels although only one type can be used at a time. The options include:

CRF-1 H00		Plate for 16 x 0.5ml CBS high security straws
CRF-1 H01		Plate for 18 x 0.3ml IMV straws
CRF-1 H02		Plate for 55 x 1.8ml cryovials (0.5ml max fill)
CRF-1 H03		Flat plate for various items/vessels
CRF-1 H05		Plate for 1 x SBS microplate
CRF-1 H06		Plate for 55 x 1.8ml cryovials (1.0ml max fill)
CRF-1 H07		Plate for 10 x 0.5ml CBS high security straws & 12 x 1.8ml cryovials (0.5ml max fill)
CRF-1 H08		Plate for 12 x 0.3ml IMV straws & 12 x 1.8ml cryovials (0.5ml max fill)



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922) 49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Ижевск (3412)26-03-58
Иваново (4932)77-34-06
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижегород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сыктывкар (8212)25-95-17
Сургут (3462)77-98-35
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольяти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Саранск (8342)22-96-24
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://grant.nt-rt.ru> || gtq@nt-rt.ru